

System Requirements

ZOLL Data Retriever for PC, version 5.10

ZOLL Data Retriever Workstation	
Computer	Pentium 2 or higher 256 MB RAM 50 MB free hard disk space
Operating system	Windows XP Professional (see Notes below)
Prerequisite	Microsoft .NET Framework 2.0
Network	Optional: Central shared folder on network to store event data files (.zol)
Peripherals (see Tested Components below)	For ZOLL AED Plus [®] and ZOLL AED Pro [®] data upload: Infrared PC port or adapter -or- USB flash drive (AED Pro only) For ZOLL M Series [®] , ZOLL E Series [®] or 1600/1700 data upload: PC card reader -or- Serial (RS-232) cable -or- Bluetooth radio (M Series and E Series only) Optional: Printer

Tested Components	
Infrared adapters	Most laptops with integrated infrared port USB infrared adapter, available from ZOLL, P/N 8000-0815
Card Reader	Most laptops with an integrated Type 2 PC card slot OmniDrive Pro external card reader (parallel), available from ZOLL, P/N 8000-0692-01 OmniDrive USB LF external card reader, available from ZOLL, P/N 8000-0695-01 OmniDrive USB2 Intern LF, internal card reader for desktop computers.
Bluetooth stack	Microsoft, WIDCOMM/BroadCOM
Serial (RS-232) cable	Serial (RS-232) cable, available from ZOLL, P/N 8000-0605-01 Ositech (RS-232 to USB) cable, available from Ositech, 1-800-563-2386 ext. 400, P/N ZOL14572-ME04.

Notes:

1. ZOLL Data Retriever is designed to function as part of the RescueNet™ Code Review system.
2. Windows 2000 Server operating system and Terminal Services (on any operating system) are **NOT** supported.
3. Any deviation from the recommended list of hardware and software cannot be officially supported. Please contact ZOLL Tech Services at (800) 348-9011 with any questions.