

Release Notes

RescueNet Code Review, ZOLL Data Relay, ZOLL Data Retriever, and Software Development Kit, Version 5.21

Document updated: October 12, 2011

© 2000 - 2011 by ZOLL Data Systems. All rights reserved. ZOLL Data Systems is a subsidiary of ZOLL Medical Corporation.

ZOLL, M Series, and E Series are registered trademarks of ZOLL Medical Corporation. RescueNet is a trademark of ZOLL Medical Corporation.

Bluetooth® is a registered trademark of Bluetooth SIG, Inc.

Other product and company names may be the trademarks of their respective owners.

Code Review

New Features

- Added new percentage to the CPR Analysis tab: % of compressions in target added to CPR Quality graph.
- For a new install of Code Review, the default CPR depth range has been updated to the 2010 Guidelines of 2 to 4 inches.
- Increased granularity of CPR rate choices in the picklist.

Bug Fixes

- Audio playback restarts the case when it reaches the end
- Updated Japanese version of the Getting Started Guide with translated screenshots.
- Updated AED Prompt translations for Japanese to match the translations spoken on the defibrillator.

Data Retriever

Bug Fixes

- Fixed time syncing when uploading via serial cable/Bluetooth.
- ZOLL Data Retriever Service: Not respecting monitoring interval configuration