

System Requirements

RescueNet™ Code Review - Standard and Enterprise Edition, version 5.10

RescueNet Code Review Workstation	
Computer	Pentium 2 or higher 256 MB RAM 50 MB free hard disk space
Operating system	Windows XP Professional (see Notes below)
Prerequisite	Microsoft .NET Framework 2.0
Network	Windows Workstation Service required to meet HIPAA standards Optional: Central shared folder on network to store event data files (.zol) Optional: Client for Microsoft Networks required if using network groups
Peripherals (see Tested Components below)	For ZOLL AED Plus® and ZOLL AED Pro® data upload: Infrared PC port or adapter -or- USB flash drive (AED Pro only) For ZOLL M Series®, ZOLL E Series® or 1600/1700 data upload: PC card reader -or- Serial (RS-232) cable -or- Bluetooth radio (M Series and E Series only) For ZOLL R Series® data upload: Compact flash card reader Optional: Printer

Tested Components	
Infrared adapters	Most laptops with integrated infrared port USB infrared adapter, available from ZOLL, P/N 8000-0815
Card Reader	Most laptops with an integrated Type 2 PC card slot OmniDrive Pro external card reader (parallel), available from ZOLL, P/N 8000-0692-01 OmniDrive USB LF external card reader, available from ZOLL, P/N 8000-0695-01 OmniDrive USB2 Intern LF, internal card reader for desktop computers.
Bluetooth stack	Microsoft, WIDCOMM/BroadCOM
Serial (RS-232) cable	Serial (RS-232) cable, available from ZOLL, P/N 8000-0605-01 Ositech (RS-232 to USB) cable, available from Ositech, 1-800-563-2386 ext. 400, P/N ZOL14572-ME04.

Notes:

1. ZOLL does not require a server class machine (hardware) for the RescueNet Code Review system.
2. Windows 2000 Server operating system and Terminal Services (on any operating system) are **NOT** supported.
3. Any deviation from the recommended list of hardware and software cannot be officially supported. Please contact ZOLL Tech Services at (800) 348-9011 with any questions.